


osom[®] COVID-19 Antigen Home Test

osom[®] COVID-19 Antigen Home Test

osom[®] COVID-19 Antigen Home Test

2 TESTS

2 TESTS



- KIT CONTENTS:**
- 2 x Individually pouched Test Devices
 - 2 x Tubes containing 300 µl of buffer liquid
 - 2 x Sterile Swabs
 - 4 x Tube holders (backside of carton)
 - 1 x Quick Reference Instructions


Rapid Results in 15 Minutes

MADE IN THE USA

CONTENTS



Single-use test device in a foil pouch



Collection tube with buffer liquid



Nasal swab package



Quick Reference Instructions



IVD



REF 1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

1067-2

30°C

86°F

15°C

59°F

15°C

59°F

15°C

59°F

15°C

59°F

15°C

59°F

15°C

59°F

15°C

59°F

15°C

59°F

15°C

59°F

15°C

59°F

Distributed by:
SEKISUI Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121 USA

3459-3

7 42860 00000 3

2 TESTS

osom[®] COVID-19 Antigen Home Test

For an Emergency Use Authorization (EUA) only

The OSOM COVID-19 Antigen Home Test is intended for the qualitative detection of SARS-CoV-2 nucleocapsid protein antigen in anterior nasal (nares) swab samples. This test does NOT determine if you had COVID-19 in the past or if you have immunity.

Determining a negative result requires multiple tests. You may need to purchase additional tests to perform serial (repeat) testing. This test is more likely to give you a false negative result when you have COVID-19 than a lab-based molecular test.

If you test positive for COVID-19 please contact your doctor/primary care physician or your local health authority immediately and adhere to the local guidelines regarding self-isolation.

For home testing (ages 2 and up).



Scan QR Code to view video instructions.

TUBE HOLDERS

A

B

C

D



LOT 42823525
(01)00742860000003
(17)241031
(10)42823525

2024-10-31

Manufactured by:
ANP Technologies, Inc.
824 Interchange Blvd.
Newark, DE 19711, USA



Item necessary to use the test but not provided in the test kit is a timer. Use within 30 minutes after opening the foil pouch. Avoid contact of the extraction liquid in Tube with skin and eyes.

SEKISUI
DIAGNOSTICS

Because every result matters™

sekisuidiagnostics.com

6 525