

OSOM - 2-PACK - JOB FORM

KIT LOT:

42825002C

EXP-26-06-30

QTY:

6400

**DON'T KIT
PRIOR TO:**

07-30-25

**CONFIRM
LABELS**



☐ _____

☐ _____

osom® COVID-19 Antigen Home 2-Test Kit

REF 1067-2 QTY 6

For Emergency Use Authorization Only



15°C 30°C
59°F 86°F

LOT 42825002C



2026-06-30

Manufactured by:
ANP Technologies, Inc.
824 Interchange Blvd.
Newark, DE 19711, USA

Distributed by:
SEKISUI Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121 USA

1907-0

☐ _____

osom® COVID-19 Antigen Home 2-Test Kit

REF 1067-2 QTY 60

For Emergency Use Authorization Only



15°C 30°C
59°F 86°F

LOT 42825002C



2026-06-30

Manufactured by:
ANP Technologies, Inc.
824 Interchange Blvd.
Newark, DE 19711, USA

Distributed by:
SEKISUI Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121 USA

CL1893-1

☐ _____

osom® Antigen Buffer Tubes

Open When Ready To Test



15°C 30°C
59°F 86°F

For use under Emergency Use Authorization (EUA) only.

LOT 252017



2026-06-30

Manufactured by:
ANP Technologies, Inc.
824 Interchange Blvd.
Newark, DE 19711, USA

Distributed by:
SEKISUI Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121 USA

1888-2

☐ _____

osom® COVID-19 Antigen Home Test Device



Open When Ready To Test

15°C 30°C

For use under Emergency Use Authorization (EUA) only.

LOT 252002



2027-01-31

Manufactured by:
ANP Technologies, Inc.
Newark, DE 19711, USA

Distributed by:
SEKISUI Diagnostics, LLC
San Diego, CA 92121 USA

☐ _____

SWAB LOT NUMBER: P1122



IMAGEN, Inc.
421 E Route 66
Glendora, CA 91740
www.imagenagency.com
Phone: 626.593.9618
Email: info@imagenagency.com

Certificate of Conformance

Product Name:	OSOM[®] COVID-19 Antigen Home Test Kit		
Product Number:	1067 - 2	Lot Number:	42825002C
Manufacture Date:		Expiration Date:	2026-06-30
Storage Condition:	15-30°C	Lot Size:	
Test Device Lot	252002	Buffer Tube Lot	252017
COC #:	COC-42825002C-0		

Test	Test Method	Acceptance Criteria	Result
Description of Contents	Visual	Kit contents include: 2 x Individually pouched Test Devices 2 x Tubes containing 300 µl of buffer liquid 2 x Sterile Swabs 4 x Tube holders (backside of carton) 1 x Quick Reference Instructions	Conforms

Inspected by: _____

Date: _____

Reviewed by: _____

Date: _____

Approved by: _____

Date: _____