

OSOM - 2-PACK - JOB FORM

KIT LOT:

42823108B

EXP-24-06-30

QTY:

11000

**DON'T KIT
PRIOR TO:**

07-30-23

IMAGEN
FOR INTERNAL USE ONLY

**CONFIRM
LABELS**



☐ _____

☐ _____

osom® COVID-19 Antigen
Home 2-Test Kit

REF 1067-2 QTY 6

For Emergency Use Authorization Only



15°C
59°F

30°C
86°F

LOT 42823108B



2024-06-30

Manufactured by:
ANP Technologies, Inc.
824 Interchange Blvd.
Newark, DE 19711, USA

Distributed by:
SEKISUI Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121 USA

1907-0

☐ _____

osom® COVID-19 Antigen
Home 2-Test Kit

REF 1067-2 QTY 60

For Emergency Use Authorization Only



15°C
59°F

30°C
86°F

LOT 42823108B



2024-06-30

Manufactured by:
ANP Technologies, Inc.
824 Interchange Blvd.
Newark, DE 19711, USA

Distributed by:
SEKISUI Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121 USA

CL1893-1

☐ _____

osom® Antigen
Buffer Tubes



Open When Ready To Test



15°C
59°F

30°C
86°F

For use under Emergency Use Authorization (EUA) only.

LOT 38423027



2024-06-30

Manufactured by:
ANP Technologies, Inc.
824 Interchange Blvd.
Newark, DE 19711, USA

Distributed by:
SEKISUI Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121 USA

1888-2

☐ _____

osom® COVID-19 Antigen Home Test
Device



Open When Ready To Test

15°C
30°C

For use under Emergency Use
Authorization (EUA) only.

LOT 41023108



2025-07-31

Manufactured by:
ANP Technologies, Inc.
Newark, DE 19711, USA

Distributed by:
SEKISUI Diagnostics, LLC
San Diego, CA 92121 USA

☐ _____

SWAB LOT NUMBER: P0985



IMAGEN, Inc.
421 E Route 66
Glendora, CA 91740
www.imagenagency.com
Phone: 626.593.9618
Email: info@imagenagency.com

Certificate of Conformance

Product Name:	OSOM[®] COVID-19 Antigen Home Test Kit		
Product Number:	1067 - 2	Lot Number:	42823108B
Manufacture Date:		Expiration Date:	2024-06-30
Storage Condition:	15-30°C	Lot Size:	
Test Device Lot	41023108	Buffer Tube Lot	38423027
COC #:	COC-42823108B-0		

Test	Test Method	Acceptance Criteria	Result
Description of Contents	Visual	Kit contents include: 2 x Individually pouched Test Devices 2 x Tubes containing 300 µl of buffer liquid 2 x Sterile Swabs 4 x Tube holders (backside of carton) 1 x Quick Reference Instructions	Conforms

Inspected by: _____

Date: _____

Reviewed by: _____

Date: _____

Approved by: _____

Date: _____