

OSOM - 2-PACK - JOB FORM

KIT LOT:

42824002


EXP-24-12-13

QTY:

13020

**DON'T KIT
PRIOR TO:**

01-13-24



**CONFIRM
LABELS**

 _____ _____

OSOM® COVID-19 Antigen Home 2-Test Kit

REF 1067-2 QTY 6

For Emergency Use Authorization Only



IVD

15°C
59°F

30°C
86°F

LOT 42824002



2024-12-13

1907-0

Manufactured by:
ANP Technologies, Inc.
824 Interchange Blvd.
Newark, DE 19711, USA

Distributed by:
SEKISUI Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121 USA

OSOM® COVID-19 Antigen Home 2-Test Kit

REF 1067-2 QTY 60

For Emergency Use Authorization Only



IVD

15°C
59°F

30°C
86°F

LOT 42824002



2024-12-13

CL1893-1

Manufactured by:
ANP Technologies, Inc.
824 Interchange Blvd.
Newark, DE 19711, USA

Distributed by:
SEKISUI Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121 USA

OSOM® Antigen Buffer Tubes

Σ 2

Open When Ready To Test



IVD

15°C
59°F

30°C
86°F

For use under Emergency Use Authorization (EUA) only.

LOT 232009



2025-06-30

1888-2

Manufactured by:
ANP Technologies, Inc.
824 Interchange Blvd.
Newark, DE 19711, USA

Distributed by:
SEKISUI Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121 USA

OSOM® COVID-19 Antigen Home Test Device



IVD

15°C
59°F

30°C
86°F

For use under Emergency Use Authorization (EUA) only.

LOT 242002

20NA-NA-NA

Manufactured by:
ANP Technologies, Inc.
824 Interchange Blvd.
Newark, DE 19711, USA

Distributed by:
SEKISUI Diagnostics, LLC
6659 Top Gun Street
San Diego, CA 92121 USA

SWAB LOT NUMBER: P1118

 _____ _____



IMAGEN, Inc.
421 E Route 66
Glendora, CA 91740
www.imagenagency.com
Phone: 626.593.9618
Email: info@imagenagency.com

Certificate of Conformance

Product Name:	OSOM® COVID-19 Antigen Home Test Kit		
Product Number:	1067 - 2	Lot Number:	42824002
Manufacture Date:		Expiration Date:	2024-12-13
Storage Condition:	15-30°C	Lot Size:	
Test Device Lot	242002	Buffer Tube Lot	232009
COC #:	COC-42824002-0		

Test	Test Method	Acceptance Criteria	Result
Description of Contents	Visual	Kit contents include: 2 x Individually pouched Test Devices 2 x Tubes containing 300 µl of buffer liquid 2 x Sterile Swabs 4 x Tube holders (backside of carton) 1 x Quick Reference Instructions	Conforms

Inspected by: _____ Date: _____

Reviewed by: _____ Date: _____

Approved by: _____ Date: _____