

OSOM - 2-PACK - JOB FORM

KIT LOT:

42824009

EXP-24-12-13

QTY:

12500

**DON'T KIT
PRIOR TO:**

01-13-24

**CONFIRM
LABELS**

(01)00742860000003
(17)241213
(10)42824009

LOT 42824009

2024-12-13



OSOM[®] COVID-19 Antigen Home 2-Test Kit

REF 1067-2 **QTY** 6

For Emergency Use Authorization Only

IVD 15°C 30°C
59°F 86°F

LOT 42824009

2024-12-13

Manufactured by: ANP Technologies, Inc. 824 Interchange Blvd. Newark, DE 19711, USA

Distributed by: SEKISUI Diagnostics, LLC 6659 Top Gun Street San Diego, CA 92121 USA

1907-0

OSOM[®] COVID-19 Antigen Home 2-Test Kit

REF 1067-2 **QTY** 60

For Emergency Use Authorization Only

IVD 15°C 30°C
59°F 86°F

LOT 42824009

2024-12-13

Manufactured by: ANP Technologies, Inc. 824 Interchange Blvd. Newark, DE 19711, USA

Distributed by: SEKISUI Diagnostics, LLC 6659 Top Gun Street San Diego, CA 92121 USA

CL1893-1

OSOM[®] Antigen Buffer Tubes

Open When Ready To Test

IVD 15°C 30°C
59°F 86°F

For use under Emergency Use Authorization (EUA) only.

LOT 242005

2025-06-30

Manufactured by: ANP Technologies, Inc. 824 Interchange Blvd. Newark, DE 19711, USA

Distributed by: SEKISUI Diagnostics, LLC 6659 Top Gun Street San Diego, CA 92121 USA

1888-2

OSOM[®] COVID-19 Antigen Home Test Device

Open When Ready To Test

IVD 15°C 30°C

For use under Emergency Use Authorization (EUA) only. **LOT** 242009

2026-01-31

Manufactured by: ANP Technologies, Inc. Newark, DE 19711, USA

Distributed by: SEKISUI Diagnostics, LLC San Diego, CA 92121 USA

SWAB LOT NUMBER: P1118



IMAGEN, Inc.
421 E Route 66
Glendora, CA 91740
www.imagenagency.com
Phone: 626.593.9618
Email: info@imagenagency.com

Certificate of Conformance

Product Name:	OSOM[®] COVID-19 Antigen Home Test Kit		
Product Number:	1067 - 2	Lot Number:	42824009
Manufacture Date:		Expiration Date:	2024-12-13
Storage Condition:	15-30°C	Lot Size:	
Test Device Lot	242009	Buffer Tube Lot	242005
COC #:	COC-42824009-0		

Test	Test Method	Acceptance Criteria	Result
Description of Contents	Visual	Kit contents include: 2 x Individually pouched Test Devices 2 x Tubes containing 300 µl of buffer liquid 2 x Sterile Swabs 4 x Tube holders (backside of carton) 1 x Quick Reference Instructions	Conforms

Inspected by: _____

Date: _____

Reviewed by: _____

Date: _____

Approved by: _____

Date: _____